



## Goat anti-human IgG (H+L), HRP Conjugated

Catalog	PMK-014-098M	PMK-014-098S	Tel :027-87888039
Quantity	100μL	1ml	E-mail: postmaster@biopmk.com
			Web:www.biopmk.com

**For research use only.**

Applications	Species Cross-Reactivity	Molecular Weight	Isotype
WB, IHC	Human	N/A	IgG1

**Storage Buffer & Condition:** PBS, pH 7.4, 15mg/ml Bovine Serum Albumin as Stabilizer and 50% Glycerol. Store at -20°C. Do not aliquot the antibody.

**Suggested Dilution Range:** 1:500-1:5,000 for immunohisto/cytochemistry  
 1:3,000-1:100,00 for ELLSA and Western blotting with chromogenic substrates  
 1:5,000-1:100,00 for Western blotting with ECL substrates

**Optimal dilutions should be determined by the end user.**

**Purity:** The antibody reacts with whole molecule human IgG. It also reacts with the light chains of other human immunoglobulins.

**Antibody Specificity:** Based on immunoelectrophoresis and/or ELISA, the antibody reacts with whole molecule mouse IgG. It also reacts with the light chains of other mouse immunoglobulins. No antibody was detected against non-immunoglobulin serum proteins.

The antibody may cross-react with immunoglobulins from other species.

**Expiration date:** Store at -20°C. Avoid repeated freeze / thaw cycles. The expiration date may be extended if test results are acceptable for the intended use.

**Note:** For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.