



Goat anti-Rat IgG (H+L), HRP Conjugated

Catalog PMK-014-099M PMK-014-099S

Tel :027-87888039

Quantity 100 μ L 1ml

E-mail: postmaster@biopmk.com

Web:www.biopmk.com

For research use only.

Applications	Species Cross-Reactivity	Molecular Weight	Isotype
WB, IHC	Rat	N/A	IgG1

Storage Buffer & Condition: PBS, pH 7.4, 15mg/ml Bovine Serum Albumin as Stabilizer and 50% Glycerol. Store at -20°C. Do not aliquot the antibody.

Suggested Dilution Range: 1:500-1:5,000 for immunohisto/cytochemistry

1:3,000-1:100,00 for ELLSA and Western blotting with chromogenic

substrates 1:50,000-1:100,00 for Western blotting with ECL substates

Optimal dilutions should be determined by the end user.

Purity:The antibody was purified from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.

Antibody Specificity:Based on immunoelectrophoresis and/or ELISA,the antibody reacts with whole molecule rat IgG. It also reacts with the light chains of other rat immunoglobulins. No antibody was detected against non-immunoglobulin serum proteins.

The antibody may cross-react with immunoglobulins from other species.

Expiration date:Store at -20°C. Avoid repeated freeze / thaw cycles.The expiration date may be extended if test results are acceptable for the intended use.

Note:For in vitro research use only,not for diagnostic or therapeutic use.This product is not a medical device.