



Goat Anti-Chicken IgY (H+L) , HRP Conjugated

Catalog	PMK-015-098M	PMK-015-098S	Tel :027-87888039
Quantity	100μL	1ml	E-mail: postmaster@biopmk.com
			Web:www.biopmk.com

For research use only.

Applications	Species Cross-Reactivity	Molecular Weight	Isotype
WB、IHC	chicken	N/A	IgG

Storage Buffer & Condition: PBS, pH 7.4, 15mg/ml Bovine Serum Albumin as Stabilizer and 50% Glycerol. Store at -20°C. Do not aliquot the antibody.

Suggested Dilution Range: 1:200-1:5,000 for immunohisto/cytochemistry
 1:3,000-1:100,000 for ELLSA and Western blotting with chromogenic substrates
 1:5,000-1:10,000 for Western blotting with ECL substates

Optimal dilutions should be determined by the end user.

Purity:The antibody was purified from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.

Antibody Specificity:Based on immunoelectrophoresis and/or ELISA,the antibody reacts with whole molecule chicken IgY. It also reacts with the light chains of other chicken immunoglobulins. No antibody was detected against non-immunoglobulin serum proteins.The antibody may cross-react with immunoglobulins from other species.

Expiration date:Store at -20°C. Avoid repeated freeze / thaw cycles.The expiration date may be extended if test results are acceptable for the intended use.

Note:For in vitro research use only,not for diagnostic or therapeutic use.This product is not a medical device.