

HA-Tag Rabbit Polyclonal Antibody

Catalog	PMK103M	PMK103S	Tel :400-457-3801
Quantity	50 μ L	100 μ L	E-mail: postmaster@biopmk.com
			Web: www.biopmk.com

For research use only.

Applications	Species Cross-Reactivity	Molecular Weight	Isotype
WB	N/A	N/A	IgG


Storage Buffer & Condition: PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
 Store at -20°C. Do not aliquot the antibody.

Recommended dilutions: WB: 1:5,000

Optimal dilutions should be determined by the end user.

Specificity: The HA tag antibody can recognize C-terminal, internal, and N-terminal HA-tag fusion proteins.

Background: Human influenza hemagglutinin (HA) is a surface glycoprotein required for the infectivity of the human virus. The HA tag is derived from the HA-molecule corresponding to amino acids 98-106 has been extensively used as a general epitope tag in expression vectors. Many recombinant proteins have been engineered to express the HA tag, which does not appear to interfere with the bioactivity or the biodistribution of the recombinant protein. HA tag antibody is a highly sensitive and affinity monoclonal antibody applicable to HA-tagged fusion protein detection.



2 μ g HA fusion protein+ Primary antibody dilution at
 1、 1:5,000
 2、 1:10,000

Application: WB-Western blot IHC-Immunochemistry IF-Immunofluorescence IP-Immunoprecipitation ChIP-Chormatin Immunoprecipitation

Reactivity: H-Human R-Rat M-Mouse Mk-Monkey Dg-Dog Ch-Chicken Hm-Hamster Rb-Rabbit Sh-Sheep Pg-Pig