

Recombinant Human Wnt3a V2

Catalog No : PMK2142

Known As: MGC119418; MGC119419; MGC119420; protein Wnt-3a; wingless-type MMTV integration site family, member 3A; Wnt3a; Wnt-3a

PROPERTIES

Description	Recombinant Human Protein Wnt-3a is produced by our Mammalian expression system and the target gene encoding Ser19-Lys352 is expressed with a fusion design at the N-terminus.
Accession	P56704
Formulation	Lyophilized from a 0.2 µm filtered solution of 10mM PB, 5% Sucrose, 0.01% Tween 80, pH7.4.
Size	10µg/50µg/500µg/1mg
Purity	> 95%
Endotoxin	< 0.01EU/µg as determined by LAL test.
Predicted Mol Mass	105.7 kDa
Apparent Mol Mass	110-155 kDa, reducing conditions
Reconstitution	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100µg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.
Shipping	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature listed below.
Storage	Lyophilized protein should be stored at ≤ -20°C, stable for one year after receipt. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at ≤ -20°C for 3 months.
Background	Wnt-3a is one of 19 vertebrate members of the Wingless-type MMTV integration site (Wnt) family of highly conserved cysteine-rich secreted glycoproteins important for normal developmental processes. Required for normal embryonic mesoderm development and formation of caudal somites. Required for normal morphogenesis of the developing neural tube (By similarity). Mediates self-renewal of the stem cells at the bottom of intestinal crypts (in vitro).

NOTE: The product listed herein is for research use only and is not intended for use in human or clinical diagnosis. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license. We cannot be responsible for patent infringements or other violations that may occur with the use of this product.